

Data Sheet (15.02.2016)

mi-Taq Set

Thermostable DNA Polymerase

Source: *Thermus aquaticus*, gene expressed in *E. coli*

Cat.-No.	Size	Conc.
mi-E6003	500 units	5 units/μl

For research use only! Only for in vitro use!

Content

- mi-Taq DNA Polymerase, 500 units (5 u/μl)
- dNTP mixture, 800 μl (2.5 mM each, in water (sodium salts, pH 7-9) Purity: ≥ 98% for each dNTP)
- 10x Buffer (MgCl₂ 15 mM), 1000 μl
- 6x Loading dye, 500 μl

Unit definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmol of dNTP into an acid-insoluble form in 30 minutes at 72°C.

Taq Pol in storage buffer

20 mM Tris-HCl (pH 8.0 at 25°C), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 50% Glycerol, 0.5% Tween 20 and 0.5% NP-40

Recommended PCR Assay (50 μl volume)

5 μl	10x Reaction buffer with MgCl ₂
0.5 μl	Taq Pol (2.5 u)
4 μl	of dNTP Mix (2.5 mM each)
0.2-1 μM	of each Primer
2-50 ng	template DNA
Fill up to 50 μl	PCR grade H ₂ O

General reaction conditions.

94°C 3 to 5 min.	(initial denaturation)
94°C 30 sec.	(denaturation)
50-55°C* 30 sec.	25-30 cycles (annealing*)
72°C 1-2 min.**	(extension**)
72°C 5 min.	(final extension)
4-8°C indefinite to store	

* Annealing temperature: TA°C = T_m -5°C

** Amplification time: 1 sec per 60 bp approximately

Storage conditions: - 20 ± 5°C

Avoid repeated freeze/thaw cycles

Description

mi-Taq DNA Polymerase is a highly pure, thermostable recombinant DNA polymerase encoded by a modified gene from *Thermus aquaticus* and expressed in *E. coli*. Its recombinant nature ensures utmost purity, reproducibility and processivity. mi-Taq DNA Polymerase possesses a 5'→3' polymerization-dependent exonuclease activity and lacks 3'→5' exonuclease activity. The enzyme leaves a single 3'-dA nucleotide overhang that makes the products suitable for cloning by a TA vector system.

Application

routine PCR, TA cloning, RT-PCR, multiplex PCR, PCR-SSCP, PCR-RFLP and PCR-RAPD.

Performance and purity tests

- Physical Purity: mi-Taq DNA Polymerase is determined to be >90% pure as judged by SDS-PAGE gel with Coomassie brilliant blue staining.
- Endonuclease Assay: No nicking activity is detectable on an ethidium bromide-stained agarose gel, when 1 μg of supercoiled plasmid DNA (pUC18) is incubated with 5 units of mi-Taq DNA Polymerase for 8 hours at 45°C and 8 hours at 70°C in 1x reaction buffer.
- Exonuclease Assay: No exonuclease activity is detectable on an ethidium bromide-stained agarose gel, when 1 μg of lambda DNA and 1 μg lambda/Hind III DNA are incubated with 5 units of mi-Taq DNA Polymerase for 8 hours at 45°C and 8 hours at 70°C in 1x reaction buffer.
- Functional Assay: mi-Taq DNA Polymerase is tested for performance in the polymerase chain reaction using 2.5 units of the enzyme to amplify lambda DNA (0.5 kb, 1.0 kb and 3 kb) and a fragment of the human beta-globin gene (408 bp).