

mi-Taq only

Thermostable DNA Polymerase (DNA free tested)

Source: *Thermus aquaticus*, strain YT-1

Cat.-No.	Size	Conc.
mi-E8001S	200 units	5 units/μl
mi-E8001L	1000 units	5 units/μl

For research use only! Only for in vitro use!

Unit definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmol of dNTP into an acid-insoluble form in 30 minutes at 72 °C.

Taq Pol in storage buffer (red cap)

10 mM K-phosphate, pH 7.4, 0.1 mM EDTA, 50 % glycerol, 0.1 % Triton X-100, 0.1 % Tween 20

10x Reaction buffer complete KCl (black cap)

500 mM KCl, 100 mM TrisHCl pH 8.8, 0.1 % Tween-20, 15 mM MgCl₂

10x Reaction buffer complete (NH₄)₂SO₄ (green cap)

160 mM (NH₄)₂SO₄, 670 mM TrisHCl pH 8.8, 0.1 % Tween-20, 25 mM MgCl₂

10x Reaction buffer (NH₄)₂SO₄ without MgCl₂ (blue cap)

160 mM (NH₄)₂SO₄, 670 mM TrisHCl pH 8.8, 0.1 % Tween-20

MgCl₂ stock solution (yellow cap)

100 mM MgCl₂

Recommended MgCl₂ concentration: 1.5-6 mM

Recommended PCR Assay

50 μl PCR assay		
5 μl	10x Reaction buffer with MgCl ₂	green cap
0.2-0.5 μl (1-2.5 u)	Taq Pol	red cap
5 μl	of dNTP Mix (2 mM each)	
0.2-1 μM	of each Primer	
2-50 ng	Template DNA	
Fill up to 50 μl	PCR grade H ₂ O	

Description

The Thermostable DNA Polymerase (94 kDa) is an enzyme that replicates DNA at 72 °C. The enzyme catalyzes the polymerization of nucleotides into duplex DNA in 5'→3' direction in the presence of magnesium ions. It also possesses a 5'→3' polymerization-dependent exo-nuclease replacement activity. The enzyme is highly purified and free of nonspecific endo- or exonucleases. It produces single 3'-dA nucleotide overhangs and is recommended for use in routine PCR, primer extension, etc.

Sensitivity

Taq DNA polymerase effectively directs PCR with templates up to 2 kb in length. High sensitivity of PCR reactions with Taq DNA polymerase in the optimal conditions – in some reactions at least 6 DNA molecules are necessary for detection. Enzyme has a very good performance in single-copy gene PCR from genomic mammalian DNA. In contrast to this enzyme, Taq DNA polymerases from the variety of suppliers contain contaminating DNA. With DNA-contaminated Taq DNA polymerase you can observe false-positive PCR results in some cases.

QC

> 98% protein homogeneity in 10% SDS-PAGE
No detectable exo-/endonuclease activities
PCR amplification tests with different templates
PCR amplification tests without templates as Negative Control
Exonuclease efficiency test showing efficient 5'-3' Exonuclease activity

Store at - 20 ± 5°C (the enzyme is stable at room temperature at least for 3 days without any loss of activity), avoid frequent thawing and freezing.