

mi-Plasmid Midiprep Kit

Cat. No. mi-PMD25 & mi-PMD50

metabion's mi-Plasmid Midiprep Kit allows the isolation of ultrapure plasmid DNA from up to 50ml culture. Plasmid DNA purified from **metabion's** proprietary anion-exchange resin is suited for use in transfection, automated sequencing and enzymatic modification.

Downstream Applications

- Transfection
- Transformation
- Ligation and cloning
- Sequencing
- *In vitro* transcription

Product Contents

Cat. No.	mi-PMD25	mi-PMD50
Preps	25	50
VP1 Buffer	120ml	265ml
VP2 Buffer	120ml	265ml
VP3 Buffer	120ml	265ml
VPN Buffer	265ml x 2	265ml x 4
VPE Buffer	130ml	265ml
RNase A (20mg/ml)	0.600ml	1.325ml
Mini-M Column	50	100
Midi-Column	25	50

Storage Conditions

metabion's mi-Plasmid Midiprep Kit can be stored at room temperature up to 12 months. If precipitate forms in any buffer, incubate at 37°C for 30 minutes to resuspend. RNase A is preferentially stored at 4°C.

Protocol

Please read the following notes before starting the procedures.

Important Notes

- All buffers need to be mixed well before use!
- Spin RNase A solution tube before use; add all of RNase A solution into VP1 Buffer bottle and mix well; store at 4°C.
- If precipitation has formed in VP2 Buffer, incubate at 55°C for 10min (alternatively 30min at 37°C) to resuspend. Do not shake VP2 Buffer, SDS present will lead to serious foaming.
- Sit VP3 Buffer on ice before use.
- The volume of buffers VP1-3 used in the protocol is developed for 50-ml sample culture. If starting sample culture is larger than 50ml, please increase the volume of buffers VP1-3 proportionally.

*Don't risk your experiment ...
trust *metabion**

www.metabion.com

Protocol

1. Culture plasmid-containing bacterial cells in 25-50ml (high-copy-number plasmids) and 100-250ml (low-copy-number plasmids) of LB medium, respectively. Grow 12-16 hours with vigorous shaking at 37°C.
2. Harvest the bacterial cells by centrifugation at 6,000 x g for 15 minutes.
3. Equilibrate Midi-Columns by applying 3ml of 98% ethanol. Allow the column to empty by gravity flow and discard the filtrate.
4. Apply 5ml of VPN Buffer to the Midi-Column, allow it to flow through by gravity and discard the filtrate.
5. Resuspend the cell pellet in 4ml of VP1 Buffer (make sure RNase A has been added!).
The bacterial cells should be completely resuspended before adding VP2 Buffer.
6. Add 4ml of VP2 Buffer, mix gently by inverting the tube and let sit for 5 minutes.
Do not vortex; vortexing will shear genomic DNA. The lysate should be clear and viscous.
7. Add 4ml of ice-cold VP3 Buffer, mix gently by rotating.
After adding VP3 Buffer, white precipitate should form.
8. Centrifuge at 20,000 x g for 15 minutes at 4°C.
20,000 x g corresponds to 12,000 and 13,000 rpm in Beckman JA-17 and Sorvall SS-34 rotors, respectively.
9. Apply the supernatant to the Midi-Column, allow it to flow through by gravity and discard the filtrate.
10. Wash the column once with 15ml of VPN Buffer by gravity flow and discard the filtrate.
11. Apply 5ml of VPE Buffer to elute DNA by gravity flow.
12. Precipitate DNA by adding 3.75ml (0.75 volumes) of room temperature isopropanol to the eluate. Mix and centrifuge at 15,000 x g for 30 minutes at 4°C. Carefully remove the supernatant.
13. Wash the DNA pellet with 5ml of room temperature 70% ethanol and centrifuge at 15,000 x g for 10 minutes. Carefully remove the supernatant.
14. Air-dry the DNA pellet for 10 minutes and dissolve the DNA in 100µl or a suitable volume of TE or ddH₂O.
15. Some insoluble material may also elute from the column at step 11. To eliminate the insoluble material, load the dissolved DNA sample onto a Mini-M™ Column, sitting in a 1.5-ml tube, and spin at full speed in a microcentrifuge for 20 seconds. Collect the eluted DNA sample in the 1.5-ml tube.
16. Short-term (up to a few weeks), store at 4°C; long-term, store at -20°C. To avoid multiple freeze-thaw cycles, preparation of aliquots is recommended.