

Reporters

Emission colour	Fluorescent Dyes	Abs max (nm)	Em max (nm)
	Atto 390	390	476
	Atto 425	439	485
	Atto 465	453	506
	6-Fam	495	520
	Atto 495	498	526
	Atto 488	500	520
	Tet	521	542
	Atto 520	517	538
	Atto 430LS	436	545
	Joe	520	545
	Yakima Yellow®	530	550
	Atto 532	532	552
	Hex	533	559
	Atto Rho6G	533	557
	Atto 550	554	576
	Tamra	543	575
	Cy3	552	565
	DY510XL	509	590
	Atto 565	564	590
	Rox	578	604
	Texas Red®	590	615
	Atto 590	593	622
	Atto Rho13	603	627
	DY480XL	500	630
	Atto Rho14	626	646
	Atto 633	630	651
	Atto 647N	646	664
	Atto 647	647	667
	Cy5	649	670
	Atto 665	662	680
	Atto 490LS	495	658
	Atto 655	663	680
	Cy5.5	683	703

Quenchers

